

**Determination of pyrrolizidine alkaloids (PA)
in plant material by SPE-LC-MS/MS**

Method Protocol

BfR-PA-Tea-2.0/2014

CONTENT

1	Scope.....	1
2	Principle	1
3	Reagents.....	1
3.1	General.....	1
3.2	Chemicals	2
3.3	Solutions	3
3.3.1	Extraction solution.....	3
3.3.2	Aqueous ammoniacal solution for neutralisation.....	3
3.3.3	2.5 % ammonia in methanol for SPE elution (black tea and green tea)	3
3.3.4	HPLC mobile phase	3
3.3.5	Standard solution for calibration	3
4	Apparatus.....	4
5	Procedure.....	6
5.1	Sample preparation (grinding of plant material).....	6
5.2	Extraction	6
5.3	SPE-procedure	7
5.4	Reconstitution of the sample	7
6	HPLC-MS/MS analysis	7
6.1	Liquid chromatographic separation	7
6.2	Mass spectrometric operation conditions	7
6.3	Measurement	8
7	Calculation	8
7.1	Calibration function	8
7.2	Quantification	8
7.3	Reporting of results	9
8	Annex.....	10
8.1	LC-MS/MS measurement.....	10
8.2	Typical Chromatogram.....	13
8.3	Provider of PA-Standards.....	14
8.4	Flow chart of the sample preparation procedure	17

1 Scope

Pyrrolizidine alkaloids (PA) are secondary plant metabolites with carcinogenic and genotoxic properties. Currently, more than 600 PA are known. They occur in plants of families of Boraginaceae, Asteraceae and Fabaceae. The worldwide spread of these plants may lead to a contamination of herbal foodstuff, herbal medicines and animal feed (EFSA Panel on Contaminants in the Food Chain (CONTAM) 2011).

This method describes the determination of the following PA in plant material: echimidine (Em), echimidine-N-oxide (EmN), erucifoline (Er), erucifoline-N-oxide (ErN), europine (Eu), europine-N-oxide (EuN), heliotrine (Hn), heliotrine-N-oxide (HnN), intermedine (Im), intermedine-N-oxide (ImN), jacobine (Jb), jacobine-N-oxide (JbN), lasiocarpine (Lc), lasiocarpine-N-oxide (LcN), lycopsamine (La), lycopsamine-N-oxide (LaN), monocrotaline (Mc), monocrotaline-N-oxide (McN), retrorsine (Re), retrorsine-N-oxide (ReN), senecionine (Sc), senecionine-N-oxide (ScN), seneciphylline (Sp), seneciphylline-N-oxide (SpN), senecivernine (Sv), senecivernine-N-oxide (SvN), senkirkine (Sk), trichodesmine (Td).

The limits of determination and quantification for toxins determined during in-house validation are listed in the Annex 8.1.

2 Principle

A test portion of plant material is sonicated twofold in aqueous sulphuric acid solution for PA extraction. After centrifugation an aliquot of the supernatant is purified by solid-phase extraction (SPE) using reversed phase C18 material. PAs are released from the cartridge using methanol. Subsequently, the eluate is evaporated to dryness and reconstituted in methanol/water (initial HPLC conditions).

For chromatographic separation, an RP-HPLC column is used with a binary gradient. Analytes are detected by triple stage quadrupole mass spectrometry. Quantification of pyrrolizidine alkaloids is accomplished by means of matrix matched calibration.

3 Reagents

3.1 General

Please note: Since the use of this method involves reagents harmful to health, appropriate precautionary and protective measures such as avoiding skin contact and using an extractor hood must be taken.

If not specified otherwise, reagents of analytical grade and solvents suitable for HPLC-MS/MS must be used. Water must be distilled in glass vessels or demineralised before use, or must be of equivalent purity.

3.2 Chemicals

- 3.2.1 echimidine (Em)
- 3.2.2 echimidine-N-oxide (EmN)
- 3.2.3 erucifoline (Er)
- 3.2.4 erucifoline-N-oxide (ErN)
- 3.2.5 europine (Eu)
- 3.2.6 europine-N-oxide (EuN)
- 3.2.7 heliotrine (Hn)
- 3.2.8 heliotrine-N-oxide (HnN)
- 3.2.9 intermedine (Im)
- 3.2.10 intermedine-N-oxide (ImN)
- 3.2.11 jacobine (Jb)
- 3.2.12 jacobine-N-oxide (JbN)
- 3.2.13 lasiocarpine (Lc)
- 3.2.14 lasiocarpine-N-oxide (LcN)
- 3.2.15 lycopsamine (La)
- 3.2.16 lycopsamine-N-oxide (LaN)
- 3.2.17 monocrotaline (Mc)
- 3.2.18 monocrotaline-N-oxide (McN)
- 3.2.19 retrorsine (Re)
- 3.2.20 retrorsine-N-oxide (ReN)
- 3.2.21 senecionine (Sc)
- 3.2.22 senecionine-N-oxide (ScN)
- 3.2.23 seneciophylline (Sp)
- 3.2.24 seneciophylline-N-oxide (SpN)
- 3.2.25 senecivernine (Sv)
- 3.2.26 senecivernine-N-oxide (SvN)
- 3.2.27 senkirkine (Sk)
- 3.2.28 trichodesmine (Td)
- 3.2.29 formic acid 98 – 100%, e.g. Sigma-Aldrich
- 3.2.30 methanol (MeOH) in LC-MS quality, e.g. Merck LiChrosolv®
- 3.2.31 sulphuric acid 98%, e.g. Merck

- 3.2.32 ammonia 32%, e.g. Merck
- 3.2.33 ammonium formate in LC-MS quality, e.g. Fluka
- 3.2.34 acetonitrile, e.g. Merck LiChrosolv®

3.3 Solutions

3.3.1 Extraction solution

2.665 mL of sulphuric acid (H₂SO₄) (3.2.31) are filled up to 1000 mL with water. The final concentration is 0.05 M.

3.3.2 Aqueous ammoniacal solution for neutralisation

To prepare the ammoniacal solution for neutralisation of sample extracts before SPE, 5 mL of ammonia (3.2.32) are filled up to 25 mL with water.

3.3.3 2.5 % ammonia in methanol for SPE elution (SPE elution for black tea and green tea)

To prepare the 2.5 % ammonia solution in methanol 7.8 mL of ammonia (3.2.32) are filled up to 100 mL with methanol (3.2.30). The solution has to be freshly prepared per working day.

3.3.4 HPLC mobile phase

Eluent A:

315 mg ammonium formate (3.2.33) are dissolved in 5 mL of water, 1 mL of formic acid (3.2.29) is added and filled up to 1000 mL with water.

Eluent B:

315 mg ammonium formate (3.2.33) are dissolved in 5 mL of water, 1 mL of formic acid (3.2.29) is added and filled up to 1000 mL with methanol (3.2.30).

3.3.5 Standard solution for calibration

Stock solution (0.1 mg/mL):

To create a stock solution, 1 mg of a pyrrolizidine alkaloid standards are weighed using analytical balance (4.4) and filled up with acetonitrile (3.2.34) in a volumetric flask to make 10 mL. The concentration of the stock solution is 0.1 mg/mL.

Standard working solution (PA mixture, 1 µg/mL)

For preparation of the standard working solution, respective volumes of each PA stock solution (0.1 mg/mL) are pipetted into a volumetric flask and is filled up with acetonitrile (3.2.34), to obtain a concentration of 1 µg/mL.

Preparation of matrix matched standards (MMS)

For a correction of matrix effects a matrix matched calibration is used. In order to obtain the same matrix effect for MMS and for the samples the blank plant material has to be processed as described in section 5. Afterwards, MMS levels have to be prepared according to table 1.

Please note: In order to prepare of a sufficient amount of blank-extract, two times 10 mL of neutralised blank-extract (5.2) from one blank sample can be purified using two SPE-cartridges. Thereby, 2 mL of reconstituted blank-extract for the MMS preparation can be received without extraction of an additional blank sample.

table 1: Matrix matched standards

	Final PA mass concentration in calibration solution ng/mL	Final PA mass concentration µg/kg	Aliquot taken from	Aliquot Volume µL	Aliquot taken from blank plant material extract µL
MMS_1	5,0	10,0	MMS_5	20	280
MMS_2	10,0	20,0	MMS_8	20	280
MMS_3	25,0	50,0	MMS_8	20	100
MMS_4	50,0	100,0	PA-Mix	10	100
MMS_5	75,0	150,0	PA-Mix	15	185
MMS_6	100,0	200,0	PA-Mix	20	180
MMS_7	125,0	250,0	PA-Mix	25	175
MMS_8	150,0	300,0	PA-Mix	60	340

4 Apparatus

4.1 General

Usual laboratory glassware and equipment should be used and, in particular, the following:

- 4.2 **centrifugal mill** with 0.5 mm sieve, e.g. Retsch
- 4.3 **various piston pipettes and multiple dispensers**, e.g. Brand
- 4.4 **analytical balance**, capable of weighing to 0,0001 g
- 4.5 **centrifuge** for 50 mL centrifuge tubes, capable of at least 5 000 x g
- 4.6 **ultrasonic bath**
- 4.7 **overhead shaker**, e.g. Heidolph
- 4.8 **laboratory shaker**, e.g. Vortex
- 4.9 **evaporation station**, e.g. TurboVap
- 4.10 **centrifuge tube** 50 mL
- 4.11 **test tubes** 15 mL
- 4.12 **volumetric flasks**, 10 and 20 mL
- 4.13 **folded filters**, e.g. Munktell
- 4.14 **SPE cartridges**: DSC-C18 SPE (Supelco), 500 mg sorbent material, 6mL
- 4.15 **SPE vacuum chamber**
- 4.16 **Membrane filter** 0.2 μm , e.g. VWR 0,5 mL centrifugal filters, modified nylon membrane
- 4.17 **HPLC vials** 2 mL
- 4.18 **Glass inserts**, 250 μL conic for HPLC vials
- 4.19 **chromatographic column**, e.g. Thermo, Hypersil Gold®; 150 x 2.1 mm; 1,9 μm
- 4.20 **LC-MS/MS system**

5 Procedure

5.1 Sample preparation (grinding of plant material)

To determine the PA-content which is representative for the entire sample, the plant material should meet the following characteristics: uniform particle size and a homogenous distribution of PA or PA-containing material, respectively. Therefore, the entire sample material is mixed with dry ice (ratio 2:1), ground to a particle size of 0.5 mm (4.2) and homogenized for example by shaking over head (4.7). As an alternative to grinding with dry ice (excellent grinding results due to shear forces and porosity of the frozen sample material), the sample material may also be ground to a particle size of 0.25 mm if there is no considerably generation of heat.

If the test material can neither be ground with dry ice nor ground to a particle size of 0.25 mm, it is also possible to increase the weighed sample amount in contrast to 5.2 to at least 10 g (particle size 0.5 mm). In order to keep a constant ratio of sample amount to extraction volume, the used volume of extraction solution (3.3.1) needs to be increased by the same factor as the weighed sample amount.

5.2 Extraction

For the extraction of PA 2.0 g \pm 0.1 g of plant are weighed into a centrifuge tube (4.10).

- | | |
|-------------------|--|
| Extraction step 1 | For the first extraction step 20 mL of the extraction solution (3.3.1) are added to the sample. The sample material has to be wetted completely before extraction in an ultrasonic bath (4.6) for 15 min at ambient temperature. |
| Centrifugation | The sample is centrifuged for 10 min \pm 2 min at 3800 x g (4.5). The supernatant (extract 1) is transferred into a clean test tube. The sediment is used for the second extraction step. |
| Extraction step 2 | Next, 20 mL of extraction solution (3.3.1) are added to the already extracted sample. The centrifuge tube is shaken vigorously to distribute the sample (the sample can also be stirred if necessary). The sample is again extracted in the ultrasonic bath for 15 min at ambient temperature. |
| Centrifugation | The sample is centrifuged applying the above mentioned conditions. The supernatant is added to the first extract. |
| Neutralisation | The combined extracts are set to pH 7 using the neutralisation solution (3.3.2). Control of the pH value is accomplished using indicator strips. Usually, about 500 μ L to 1000 μ L of the solution are needed. |

The complete neutralised extract is passed through a folded filter (4.13). An aliquot of the filtrate is used for SPE. The filtration step prior to SPE can be repeated in case of larger quantities of remaining particles in the solution. Thereby, blockage of SPE cartridges can be avoided.

5.3 SPE-procedure

The Solid Phase Extraction (4.14) is carried out using a vacuum chamber (4.15).

Conditioning step 1	5 mL of methanol (3.2.30)
Conditioning step 2	5 mL of water
Sample load	10 mL sample (filtrate of neutralised extract)
Washing step	2 x 5 mL of water
Drying of cartridges	5 - 10 min (use the vacuum chamber (4.15))
Elution of PA	2 x 5 mL methanol (3.2.30) or in case of green and black tea 2 x 5 mL 2.5 % ammonia in methanol (3.3.3)

The eluate is dried under a nitrogen stream at $50\text{ °C} \pm 5\text{ °C}$.

5.4 Reconstitution of the sample

The residue is dissolved in 1 mL of methanol/water (5/95, v/v) by shaking (4.8).

The reconstituted sample extracts are filtered through 0.2 μm membrane filters (4.16). When using centrifugal filters, 500 μL of the sample are centrifuged at 20 000 x g for 10 min \pm 3 min. 200 μL of the filtrate are transferred into an HPLC vial (4.17) with a glass insert (4.18).

6 HPLC-MS/MS analysis

6.1 Liquid chromatographic separation

The measurements can be carried out with different high-performance liquid chromatographs (HPLC) and separation columns. The chromatographic conditions can be chosen freely. The acceptable minimum retention time is twice the retention time for the dead volume of the column. Analytes which cannot be distinguished by means of mass spectrometry must be separated chromatographically. The conditions listed in the annex 8.1 using a C18-column (4.19) and the mobile phase described in 3.3.4 have shown to be suitable in pre-trials. However, they are to be seen as examples only.

6.2 Mass spectrometric operation conditions

The measurements can be carried out with MS/MS devices of different manufacturers. In the annex 8.1, the device-specific settings of one measuring system are given as an example. These conditions have shown to be suitable in pre-trials.

Please note: For the qualitative detection and for quantification, it is necessary to detect and report two substance-specific transitions per analyte.

6.3 Measurement

For a quantitative analysis, the following criteria are defined.

Injection:

Samples and standards are injected in duplicate in order to assess repeatability of MS detection and to check for possible response drift during the sequence.

Sequence

To determine pyrrolizidine alkaloids, the following array of analysis is defined in a sequence.

1. Matrix matched standards (5 – 150 ng/mL)
2. Solvent blank
3. Samples (first injection)
4. Solvent blank
5. Matrix matched standards (5 – 150 ng/mL)
6. Solvent blank
7. Samples (second injection)

7 Calculation

The quantitative determination is performed according to the method of the matrix matched standard by integration of the peak areas in relation to the calibration line.

7.1 Calibration function

Equation1: Calibration function

$$f_{(x)} = y = ax + b$$

where

- y is the peak area of the target analyte
 a is the slope of the calibration function
 x is the concentration of the target analyte [ng/mL] in the MMS
 b is the intercept of the calibration function

7.2 Quantification

Equation2: Calculation of the PA content (analysis equation)

$$PA \text{ concentration} = \beta \times DF = \left[(y - b) \times \frac{1}{a} \right] \times \frac{V_{Extract}}{m_{weight}} \times \frac{1}{V_{Application}} \times V_{sample}$$

where

β	is the analyte concentration [ng/mL] in the sample extract
DF	is the conversion factor from ng/mL to $\mu\text{g/kg}$
y	is the peak area of the target analyte
a	is the axis intercept from the matrix calibration
a	is the increase from the matrix calibration
V_{Extract}	is the volume of extraction agent [mL]
$m_{\text{Weight-in quantity}}$	is the sample weight in [g]
$V_{\text{Application}}$	is the volume of the extract applied for SPE [mL]
V_{Sample}	is the final sample volume [mL]

7.3 Reporting of results

The results are reported in $\mu\text{g/kg}$ with two significant decimals. To convert the concentration from ng/mL injected solution to $\mu\text{g/kg}$ plant material a factor of 2 is used according to the sample preparation procedure described in chapter 5.

Reference list

DIN ISO 32645. (1994) Chemical Analysis; Decision limit, Detection limit and determination limit, Estimation in case of repeatability, terms, methods, evaluation. Deutsches Institut für Normung DIN.

EFSA Panel on Contaminants in the Food Chain (CONTAM). (2011) Scientific Opinion on Pyrrolizidine alkaloids in food and feed. The EFSA Journal 9, 1-135

8 Annex

8.1 LC-MS/MS measurement

LC-MS/MS system consisting of

Triple quadrupole mass spectrometer (Thermo TSQ Vantage)
HPLC system HPLC pump (e. g. Thermo Accela 1250),
 Degasser
 Autosampler (e. g. CTC Analytics PAL ATS MYX)
 Column oven (e. g. MayLab MistraSwitch)

HPLC settings

Eluent A refer to (3.3.4)
Eluent B refer to (3.3.4)
Column temperature 40 °C
Flow rate 300 µL/min
Injection volume 10 µL
Column e. g. Thermo Hypersil Gold; 150 x 2.1 mm, 1.9 µm
Total runtime 15 minutes

Gradient

Time (min)	% A	% B
0.0	95	5
0.5	95	5
7.0	50	50
7.5	20	80
7.6	0	100
9.0	0	100
9.1	95	5
15.0	95	5

MS settings

Ionisation Electrospray positive (ESI +)
Ion spray voltage [V] 3500 (positive polarity)
Capillary temperature [°C] 270
Vaporiser temperature [°C] 300
Sheath gas pressure [psi] 45.0
Ion sweep gas pressure [psi] 2.0
Aux gas pressure [psi] 10

Substance-specific parameters

The analytes are detected by Selected Reaction Monitoring (SRM). For analyte identification, two PA specific transitions to two product ions are chosen. The relevant transitions and the collision energy (CE) can be found in table 2. The table also lists the retention time per analyte which apply for the HPLC settings described above.

table 2: Substance-specific parameters of the LC-MS/MS method

analyte	precursor	fragment	CE	S Lens	retention time [min]
Mc	326.2	120.3	35	130	4.25
		237.3	25	130	
McN	342.1	118.3	37	141	4.99
		137.4	29	141	
ErN	366.1	136.1	30	129	5.16
		120.1	33	129	
Jb	352.1	120.1	36	110	5.25
		155.2	29	120	
Eu	330.1	138.1	20	89	5.34
		156.2	28	89	
Im	300.1	138.3	18	112	5.40
		156.3	28	112	
JbN	368.1	120.1	32	110	5.51
		296.1	23	110	
La	300.1	138.3	18	112	5.53
		156.3	28	112	
EuN	346.1	111.2	41	91	5.63
		172.1	31	91	
ImN	316.1	111.2	37	95	5.91
		138.1	26	95	
LaN	316.1	111.2	37	95	6.02
		138.1	26	95	
Td	354.2	120.3	35	137	6.37
		222.3	28	137	
ReN	368.2	136.2	30	145	6.41
		118.2	40	145	
Sp	334.2	120.3	26	138	6.56
		138.4	28	138	
Hn	314.2	138.3	19	119	6.72
		156.3	28	119	
Er	350.2	120.3	32	110	4.87
		138.1	30	110	
SpN	350.2	118.2	36	135	6.79
		136.3	32	135	
HnN	330.2	138.2	22	121	7.03
		172.1	27	121	
Sv	336.2	120.1	27	135	7.26
		138.1	27	135	
Sc	336.2	120.2	27	130	7.33
		138.2	29	130	
SvN	352.1	118.1	30	110	7.42
		120.1	36	110	
Re	352.2	120.3	27	140	7.54
		138.3	29	140	

analyte	precursor	fragment	CE	S Lens	retention time [min]
ScN	352.2	118.1	28	116	7.54
		136.3	27	116	
EmN	414.2	254.1	32	129	8.01
		352.1	27	129	
Em	398.2	120.3	23	139	8.02
		220.3	17	139	
Sk	366.2	150.3	24	132	8.19
		168.2	28	132	
Lc	412.2	120.2	30	139	8.99
		336.3	17	139	
LcN	428.2	136.1	29	135	9.33
		254.1	27	135	

table 3: Limits of detection (LOD) and limits of quantification (LOQ) determined during in-house validation of the described method*

analyte	LOD [$\mu\text{g}/\text{kg}$]	LOQ [$\mu\text{g}/\text{kg}$]
Mc	0,9	2,8
McN	1,7	5,4
ErN	1,2	3,8
Jb	1,3	4,0
Eu	0,7	2,1
Im	1,0	3,1
JbN	1,3	4,2
La	2,0	6,4
EuN	0,7	2,3
ImN	1,2	3,8
LaN	1,5	4,9
Td	1,0	3,1
ReN	1,4	4,6
Sp	1,3	4,0
Hn	0,5	1,7
Er	0,6	1,9
SpN	0,9	2,7
HnN	0,6	2,0
Sv	1,7	5,3
Sc	1,8	5,9
SvN	0,8	2,6
Re	0,8	2,7
ScN	0,9	2,9
EmN	1,9	6,1
Em	0,8	2,6
Sk	0,8	2,4
Lc	0,8	2,4
LcN	0,9	2,8

* LOD and LOQ were determined according to DIN EN ISO 32645 Calibration method (DIN ISO 32645 1994)

8.2 Typical Chromatogram




Figure 1: Typical chromatogram of a PA standard mixture (1 ng/mL), TIC of SRM-transitions

8.3 Provider of PA-Standards

pyrrolizidine alkaloid	mass	CAS	provider	order code
echimidine	397,47	520-68-3	Oskar Tropitzsch	7550006
			PhytoLab*	89553
			PlantaAnalytica	-
erucifoline	349,38	40158-95-0	Carl Roth	1657.1
			Oskar Tropitzsch	7550021
			PhytoLab*	83446
erucifoline-N-oxide	365,37	123864-94-8	Carl Roth	1664.1
			PhytoLab*	83434
europine-hydrochloride	365,86	570-19-4	Carl Roth	1676.1
			PhytoLab*	83237
europine-N-oxide	345,39	65582-53-8	AppliChem	A9574,0010
			Carl Roth	1687.1
			Oskar Tropitzsch	7500063
			PhytoLab*	83238
heliotrine	313,40	303-33-3	AppliChem	A9583,0020
			Carl Roth	1929.1
			Latoxan*	L6007
			Oskar Tropitzsch	7550511
			PhytoLab	80403
heliotrin-N-oxide	329,39	6209-65-0	AppliChem	A9590,0010
			Carl Roth	1944.1
			Oskar Tropitzsch*	755054
			PhytoLab	83236
indicine-hydrochloride	335,83	1195140-94-3	Carl Roth	1960.1
			Oskar Tropitzsch	7500069
			PhytoLab	83234
indicine-N-oxide	315,36	41708-76-3	AppliChem	A9593,0010
			Carl Roth	1961.1
			Oskar Tropitzsch	7500070
			PhytoLab	83235
intermedine	299,37	10285-06-0	Carl Roth	1962.1
			Oskar Tropitzsch	7501610
			PhytoLab*	82424
intermedine-N-oxide	315,36	95462-14-9	PhytoLab*	83434
lasiocarpine	411,49	303-34-4	AppliChem	A9596,0010
			Carl Roth	2090.1
			Oskar Tropitzsch*	7500019
			PhytoLab	80412
lasiocarpine-N-oxide	457,5	127-30-0	AppliChem	A9600,0010
			Carl Roth	2202.1
			Oskar Tropitzsch*	7501284
			PhytoLab	83220
lycopsamine	299,37	10285-07-1	Carl Roth	2208.1
			Oskar Tropitzsch	7501080
			PhytoLab*	89726
lycopsamine-N-oxide	315,36	95462-15-0	Oskar Tropitzsch PhytoLab*	7501358 83447

pyrrolizidine alkaloid	mass	CAS	provider	order code
monocrotaline	325,35	315-22-0	Carl Roth	3418.1
			Fluka	37024
			Sigma	C2401
			Oskar Tropitzsch	7550522
			PhytoLab*	89251
			R&D Chemicals	7351
Santa Cruz Biotechnology	sc-211921			
monocrotaline-N-oxide	341,36	35337-98-5	Carl Roth	2249.1
			Oskar Tropitzsch	7501658
			PhytoLab*	82629
retrorsine	351,40	480-54-6	AppliChem	A4922,0020
			Carl Roth	1213.1
			Fluka	37025
			Oskar Tropitzsch	7550659
			PhytoLab	89775
			Santa Cruz Biotechnolog	sc-215805
Sigma*	R0382			
retrorsine-N-oxide	367,40	15503-86-3	AppliChem	A8668,0010
			Carl Roth	6733.1
			Oskar Tropitzsch	7500347
			PhytoLab*	82630
senecionine	335,40	130-01-8	AppliChem	A2071,0020
			Carl Roth*	2261.1
			Fluka	37031
			Oskar Tropitzsch	7550292
			PhytoLab*	89789
			R&D Chemicals	1828
Sigma	17806			
Santa Cruz Biotechnology	sc-286770			
senecionine-N-oxide	351,40	13268-67-2	AppliChem	A8678,0010
			Carl Roth	6734.1
			Oskar Tropitzsch	7500301
			PhytoLab*	82631
seneciphylline	333,39	480-81-9	AppliChem	A2072,0020
			Carl Roth*	6414.1
			Fluka	37033
			R&D Chemicals	1850
			Santa Cruz Biotechnology Inc.	sc-229697
			ABCR GmbH	AB167974
PhytoLab*	89275			
seneciphylline-N-oxide	349,38	38710-26-8	AppliChem	A8684,0010
			Carl Roth	6735.1
			Oskar Tropitzsch	7500573
			PhytoLab*	82632
senecivernine	335.40	72755-25-0	Carl Roth	2209.1
			Oskar Tropitzsch	7550066
			PhytoLab*	83436

pyrrolizidine alkaloid	mass	CAS	provider	order code
senecivernine-N-oxide	351,39	101687-28-9	Carl Roth PhytoLab*	2215.1 83437
senkirkine	365,43	2318-18-5	AppliChem	A6765,0010
			Carl Roth	4934.1
			Fluka	37032
			Oskar Tropitzsch	7500441
trichodesmine	353,41	548-90-3	PhytoLab* Latoxan*	89274 L6049

* substances were used for the in-house validation by BfR

8.4 Flow chart of the sample preparation procedure

