

Quantifying the bias of the viable cell enumeration process and its impact on microbial inactivation




International Conference on Uncertainty in Risk Assessment

Alberto Garre, Jose A. Egea, Paula M. Periago, Alfredo Palop, Pablo S. Fernandez

Universidad Politécnica de Cartagena (Spain)

Serial dilution/plate count




Sources of uncertainty




- (Sub)sampling error.
- Log-transformation.
- Errors in volumes (calibration errors).
- Counting errors.
- Clustering of bacterial cells.
- ...



Modelling of one (sub)sampling




- Finite volumes.
- Finite number of particles.
- Uncorrelated positions (no clumps; perfect mixing).




$$N_1 \sim \text{Binom}(N_0 \cdot V_0, p = \frac{V_1}{V_0})$$



Modelling of serial dilutions




$$N_i | N_{i-1} \sim \text{Binom}(N_{i-1} \cdot V_{i-1}, p = \frac{V_i}{V_{i-1}})$$



$$N_n | N_0 \sim \text{Binom}(N_0 \cdot V_0, p = \frac{V_n}{V_0} = f^n)$$



Case study: microbial inactivation



Survivor curve typically non-linear

Shoulder: Se debe de mantener el tratamiento para inactivar

Tail: More resistant fraction/acclimation.

Very relevant for food safety

How are non-linearities affected by experimental error?




Simulation of inactivation experiments



1. Log-linear survivor curve $\log N(t) = \log N_0 - \frac{1}{D} t$
2. 11 sampling points. Dilutions 0 to -5.
3. Plate count according to Binomial distribution.
4. Plates with >300 colonies omitted.
5. 600 experiments simulated.



Uncertainty can cause bias



Experimental design affects bias



	0.1 mL en placa	0.5 mL en placa	1 mL en placa
0 dilutions	0.87 log CFU/ml	0.17 log CFU/ml	-0.17 log CFU/ml
1 dilution	1.86 log CFU/ml	1.17 log CFU/ml	0.84 log CFU/ml
2 dilutions	2.89 log CFU/ml	2.17 log CFU/ml	1.88 log CFU/ml
3 dilutions	3.94 log CFU/ml	3.22 log CFU/ml	2.92 log CFU/ml
4 dilutions	5.00 log CFU/ml	4.25 log CFU/ml	3.93 log CFU/ml




Conclusions/outlook



- Uncertainty related to experimental design.
- Uncertainty can cause bias.
- Reporting of volumes plated/dilutions in scientific publications.
- Do not mix different dilutions/volumes plated.
- Detection limit (Chik et al. 2018, *Learning Something From Nothing: The Critical Importance of Rethinking Microbial Non-detects*. *Frontiers in Microbiology*).





Food Research International xxx (2019) xxx-xxx



Contents lists available at ScienceDirect

Food Research International

journal homepage: www.elsevier.com



Tail or artefact? Illustration of the impact that uncertainty of the serial dilution and cell enumeration methods has on microbial inactivation


Alberto Garre^a, Jose A. Egea^b, Arturo Esnoz^a, Alfredo Palop^a, Pablo S. Fernandez^{a, *}

^a *Departamento de Ingeniería de Alimentos y del Equipamiento Agrícola, Instituto de Biotecnología Vegetal, Universidad Politécnica de Cartagena (ETSIA), Paseo Alfonso XIII, 48, 30203 Cartagena, Spain*

^b *Centro de Edafología y Biología Aplicada del Segura (CEBAS-CSIC), Campus Universitario de Espinardo, E-30100, Murcia, Spain*




Acknowledgements



- The financial support of this research work was provided by the Ministry of Economy and Competitiveness (MINECO) of the Spanish Government and European Regional Development Fund (ERDF) through projects AGL2013-48993-C2-1-R and DPI2011-28112-C04-04.
- Alberto Garre (BES-2014-070946) is grateful to the MINECO for awarding him a pre-doctoral grant



Danke schön!
¡Muchas gracias!



alberto.garre@upct.es